

## Compliance Chart

Inflation Pressure [bar]	Stent Outer Diameter [mm]					
	Ø 5.0	Ø 6.0	Ø 7.0	Ø 8.0	Ø 9.0	Ø 10.0
8				8.0	9.0	10.0
9	5.0	6.0	7.0	8.3	9.2	10.2
10	5.2	6.2	7.2	8.5	9.4	10.4
11	5.3	6.4	7.3	8.7	9.6	10.6
12	5.4	6.5	7.5	8.8	9.7	10.7
13	5.5	6.6	7.6			

**NP** Nominal Pressure

**RBP** Rated Burst Pressure