

Compliance Chart

Inflation Pressure [bar]	Stent Outer Diameter [mm]					
	Ø 5.0	Ø 6.0	Ø 7.0	Ø 8.0	Ø 9.0	Ø 10.0
8			7.0	8.0	9.0	10.0
9			7.2	8.2	9.2	10.3
10	5.0	6.0	7.4	8.4	9.4	10.5
11	5.2	6.2	7.5	8.5	9.6	10.6
12	5.3	6.3	7.6	8.6	9.8	10.8
13	5.4	6.4				

NP

Nominal Pressure

RBP

Rated Burst Pressure